

产品名称: KBTBD8 Mouse Monoclonal Antibody
产品货号: AMM81455



产品概述 (Summary)

产品名称 (Production Name)	KBTBD8 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,ELISA,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	KBTBD8
别名 (Alternative Names)	TA-KRP
基因 ID (Gene ID)	84541.0
蛋白 ID (SwissProt ID)	Q8NFY9.Purified recombinant fragment of human KBTBD8 (AA: 264-464) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	68.8kDa

研究背景 (Background)

KBTBD8 (kelch repeat and BTB (POZ) domain containing 8) is a protein-coding gene. An important paralog of this gene is


产品名称: KBTBD8 Mouse Monoclonal Antibody
产品货号: AMM81455




KLHL21.

研究领域 (Research Area)


图片 (Image Data)



Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);



Western blot analysis using KBTBD8 mouse mAb against A431 (1), Jurkat (2), Hela (3), K562 (4), and HEK293 (5) cell lysate.



Flow cytometric analysis of A431 cells using KBTBD8 mouse mAb (green) and negative control (red).

产品名称: KBTBD8 Mouse Monoclonal Antibody
产品货号: AMM81455



注意事项 (Note)

For research use only .