

产品名称: EMD Mouse Monoclonal Antibody
产品货号: AMM81475



产品概述 (Summary)

产品名称 (Production Name)	EMD Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	IHC, ICC, ELISA, FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	EMD
别名 (Alternative Names)	STA; EDMD; LEMD5
基因 ID (Gene ID)	2010.0
蛋白 ID (SwissProt ID)	P50402. Purified recombinant fragment of human EMD (AA: 1-222) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:200-1:1000, ICC 1:200-1:1000, ELISA 1:5000-1:20000, FC 1:200-1:400
蛋白分子量 (Molecular Weight)	29kDa

研究背景 (Background)

Emerin is a serine-rich nuclear membrane protein and a member of the nuclear lamina-associated protein family. It


产品名称: EMD Mouse Monoclonal Antibody
产品货号: AMM81475




mediates membrane anchorage to the cytoskeleton. Dreifuss-Emery muscular dystrophy is an X-linked inherited degenerative myopathy resulting from mutation in the emerin gene.

研究领域 (Research Area)


图片 (Image Data)



Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);




Immunofluorescence analysis of HeLa cells using EMD mouse mAb (green). Blue: DRAQ5 fluorescent DNA dye. Red: Actin filaments have been labeled with Alexa Fluor- 555 phalloidin. Secondary antibody from Fisher (Cat#: 35503)




Flow cytometric analysis of A549 cells using EMD mouse mAb (green) and negative control (red).

产品名称: EMD Mouse Monoclonal Antibody
产品货号: AMM81475



Flow cytometric analysis of K562 cells using EMD mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human ovarian cancer tissues using EMD mouse mAb with DAB staining.

注意事项 (Note)

For research use only .