

产品名称: NFKBIA Mouse Monoclonal Antibody
产品货号: AMM81483



产品概述 (Summary)

产品名称 (Production Name)	NFKBIA Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	ICC,ELISA,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	NFKBIA
别名 (Alternative Names)	IKBA; MAD-3; NFKBI
基因 ID (Gene ID)	4792.0
蛋白 ID (SwissProt ID)	P25963.Purified recombinant fragment of human NFKBIA (AA: 150-291) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	ICC 1:200-1:1000,ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	35.6kDa

研究背景 (Background)

This gene encodes a member of the NF-kappa-B inhibitor family, which contain multiple ankrin repeat domains. The


产品名称: NFKBIA Mouse Monoclonal Antibody
产品货号: AMM81483




encoded protein interacts with REL dimers to inhibit NF-kappa-B/REL complexes which are involved in inflammatory responses. The encoded protein moves between the cytoplasm and the nucleus via a nuclear localization signal and CRM1-mediated nuclear export. Mutations in this gene have been found in ectodermal dysplasia anhidrotic with T-cell immunodeficiency autosomal dominant disease.

研究领域 (Research Area)

图片 (Image Data)




Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);




Immunofluorescence analysis of A549 cells using NFKBIA mouse mAb (green). Blue: DRAQ5 fluorescent DNA dye. Red: Actin filaments have been labeled with Alexa Fluor- 555 phalloidin. Secondary antibody from Fisher (Cat#: 35503)

产品名称: NFKBIA Mouse Monoclonal Antibody
产品货号: AMM81483



Immunofluorescence analysis of HeLa cells using NFKBIA mouse mAb (green). Blue: DRAQ5 fluorescent DNA dye. Red: Actin filaments have been labeled with Alexa Fluor- 555 phalloidin. Secondary antibody from Fisher (Cat#: 35503)



Flow cytometric analysis of A549 cells using NFKBIA mouse mAb (green) and negative control (red).

注意事项 (Note)

For research use only .