

产品名称: THAP1 Mouse Monoclonal Antibody
产品货号: AMM81502



产品概述 (Summary)

产品名称 (Production Name)	THAP1 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	ICC,ELISA,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	THAP1
别名 (Alternative Names)	DYT6
基因 ID (Gene ID)	55145.0
蛋白 ID (SwissProt ID)	Q9NVV9.Purified recombinant fragment of human THAP1 (AA: 1-213) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	ICC 1:200-1:1000,ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	24.9kDa

研究背景 (Background)

The protein encoded by this gene contains a THAP domain, a conserved DNA-binding domain. This protein colocalizes with


产品名称: THAP1 Mouse Monoclonal Antibody
产品货号: AMM81502




the apoptosis response protein PAWR/PAWR-4 in promyelocytic leukemia (PML) nuclear bodies, and functions as a proapoptotic factor that links PAWR to PML nuclear bodies. Alternatively spliced transcript variants encoding distinct isoforms have been observed.

研究领域 (Research Area)


图片 (Image Data)



Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);



Immunofluorescence analysis of HeLa cells using THAP1 mouse mAb (green). Blue: DRAQ5 fluorescent DNA dye. Red: Actin filaments have been labeled with Alexa Fluor- 555 phalloidin. Secondary antibody from Fisher (Cat#: 35503)



产品名称: THAP1 Mouse Monoclonal Antibody
产品货号: AMM81502



Flow cytometric analysis of Jurkat cells using THAP1 mouse mAb (green) and negative control (red).

注意事项 (Note)

For research use only .