

产品名称: RAB11FIP1 Mouse Monoclonal Antibody
产品货号: AMM81508



产品概述 (Summary)

产品名称 (Production Name)	RAB11FIP1 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IHC,ELISA,FC
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	RAB11FIP1
别名 (Alternative Names)	RCP; NOEL1A; rab11-FIP1
基因 ID (Gene ID)	80223.0
蛋白 ID (SwissProt ID)	Q6WKZ4 .Purified recombinant fragment of human RAB11FIP1 (AA: 130-271) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:200-1:1000,ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	137.2kDa

研究背景 (Background)

This gene encodes one of the Rab11-family interacting proteins (Rab11-FIPs), which play a role in the Rab-11 mediated


产品名称: RAB11FIP1 Mouse Monoclonal Antibody
产品货号: AMM81508




recycling of vesicles. The encoded protein may be involved in endocytic sorting, trafficking of proteins including integrin subunits and epidermal growth factor receptor (EGFR), and transport between the recycling endosome and the trans-Golgi network. Alternative splicing results in multiple transcript variants. A pseudogene is described on the X chromosome

研究领域 (Research Area)


图片 (Image Data)



Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);




Western blot analysis using RAB11FIP1 mouse mAb against Raji (1), A431 (2), HCT116 (3), LOVO (4), SW620 (5), SW480 (6), HepG2 (7), and NIH3T3 (8) cell lysate.




产品名称: RAB11FIP1 Mouse Monoclonal Antibody
产品货号: AMM81508



Flow cytometric analysis of Raji cells using RAB11FIP1 mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human cervical cancer tissues using RAB11FIP1 mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human ovarian cancer tissues using RAB11FIP1 mouse mAb with DAB staining.

注意事项 (Note)

For research use only .