

产品名称: LRP1 Mouse Monoclonal Antibody
产品货号: AMM81532



产品概述 (Summary)

产品名称 (Production Name)	LRP1 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	IHC, ICC, ELISA, FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	LRP1
别名 (Alternative Names)	APR; LRP; A2MR; CD91; APOER; LRP1A; TGFBR5; IGFBP3R
基因 ID (Gene ID)	4035.0
蛋白 ID (SwissProt ID)	Q07954. Purified recombinant fragment of human LRP1 (AA: 20-155) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:200-1:1000, ICC 1:200-1:1000, ELISA 1:5000-1:20000, FC 1:200-1:400
蛋白分子量 (Molecular Weight)	50.5kDa

研究背景 (Background)

The protein encoded by this gene is an endocytic receptor involved in several cellular processes, including intracellular


产品名称: LRP1 Mouse Monoclonal Antibody
产品货号: AMM81532




signaling, lipid homeostasis, and clearance of apoptotic cells. In addition, the encoded protein is necessary for the A2M-mediated clearance of secreted amyloid precursor protein and beta-amyloid, the main component of amyloid plaques found in Alzheimer patients. Expression of this gene decreases with age and has been found to be lower than controls in brain tissue from Alzheimer patients.

研究领域 (Research Area)


图片 (Image Data)



Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);




Flow cytometric analysis of Jurkat cells using LRP1 mouse mAb (green) and negative control (red).



产品名称: LRP1 Mouse Monoclonal Antibody
产品货号: AMM81532



Immunohistochemical analysis of paraffin-embedded human esophageal cancer tissues using LRP1 mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human stomach cancer tissues using LRP1 mouse mAb with DAB staining.

注意事项 (Note)

For research use only .