

产品名称: PDGFRB Mouse Monoclonal Antibody
产品货号: AMM81534



产品概述 (Summary)

产品名称 (Production Name)	PDGFRB Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	ELISA,FC
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PDGFRB
别名 (Alternative Names)	IMF1; IBGC4; JTK12; PDGFR; CD140B; PDGFR1; PDGFR-1
基因 ID (Gene ID)	5159.0
蛋白 ID (SwissProt ID)	P09619. Purified recombinant fragment of human PDGFRB (AA: 33-133) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	124kDa

研究背景 (Background)

This gene encodes a cell surface tyrosine kinase receptor for members of the platelet-derived growth factor family. These

产品名称: PDGFRB Mouse Monoclonal Antibody
产品货号: AMM81534




growth factors are mitogens for cells of mesenchymal origin. The identity of the growth factor bound to a receptor monomer determines whether the functional receptor is a homodimer or a heterodimer, composed of both platelet-derived growth factor receptor alpha and beta polypeptides. This gene is flanked on chromosome 5 by the genes for granulocyte-macrophage colony-stimulating factor and macrophage-colony stimulating factor receptor; all three genes may be implicated in the 5-q syndrome. A translocation between chromosomes 5 and 12, that fuses this gene to that of the translocation, ETV6, leukemia gene, results in chronic myeloproliferative disorder with eosinophilia.


研究领域 (Research Area)

PI3K-Akt signaling pathway, MAPK signaling pathway, Hippo signaling pathway

图片 (Image Data)



Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);



Flow cytometric analysis of NIH/3T3 cells using PDGFRB mouse mAb (green) and negative control (red).

注意事项 (Note)

For research use only .