

产品名称: HIST2H3C(27Ac) Mouse Monoclonal Antibody
产品货号: AMM81576



产品概述 (Summary)

产品名称 (Production Name)	HIST2H3C(27Ac) Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	IHC,ELISA,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	HIST2H3C(27Ac)
别名 (Alternative Names)	H3; H3.2; H3/M; H3F2; H3FM; H3FN
基因 ID (Gene ID)	126961.0
蛋白 ID (SwissProt ID)	Q71DI3.Synthesized peptide of human HIST2H3C (AA: ATKAARK(Ac)SAPATGGV) .

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:200-1:1000,ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	15.4kDa

研究背景 (Background)

Histones are basic nuclear proteins that are responsible for the nucleosome structure of the chromosomal fiber in


产品名称: HIST2H3C(27Ac) Mouse Monoclonal Antibody
产品货号: AMM81576




eukaryotes. This structure consists of approximately 146 bp of DNA wrapped around a nucleosome, an octamer composed of pairs of each of the four core histones (H2A, H2B, H3, and H4). The chromatin fiber is further compacted through the interaction of a linker histone, H1, with the DNA between the nucleosomes to form higher order chromatin structures. This gene is intronless and encodes a member of the histone H3 family. Transcripts from this gene lack polyA tails; instead, they contain a palindromic termination element. This gene is found in a histone cluster on chromosome 1. This gene is one of four histone genes in the cluster that are duplicated; this record represents the telomeric copy.

研究领域 (Research Area)

图片 (Image Data)




Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);




Flow cytometric analysis of *** cells using *** mouse mAb (green) and negative control (red).


产品名称: HIST2H3C(27Ac) Mouse Monoclonal Antibody
产品货号: AMM81576



Flow cytometric analysis of *** cells using HIST2H3C(27Ac) mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human colon cancer tissues using HIST2H3C(27Ac) mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human breast cancer tissues using HIST2H3C(27Ac) mouse mAb with DAB staining.

注意事项 (Note)

For research use only .