

产品名称: BMP4 Mouse Monoclonal Antibody
产品货号: AMM81584



产品概述 (Summary)

产品名称 (Production Name)	BMP4 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,ICC,ELISA
种属反应性 (Reactivity)	Human,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	BMP4
别名 (Alternative Names)	ZYME; BMP2B; OFC11; BMP2B1; MCOPS6
基因 ID (Gene ID)	652.0
蛋白 ID (SwissProt ID)	P12644.Purified recombinant fragment of human BMP4 (AA: 277-408) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ICC 1:100-1:500,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	46.5kDa

研究背景 (Background)

The protein encoded by this gene is a member of the bone morphogenetic protein family which is part of the transforming

产品名称: BMP4 Mouse Monoclonal Antibody
产品货号: AMM81584




growth factor-beta superfamily. The superfamily includes large families of growth and differentiation factors. Bone morphogenetic proteins were originally identified by an ability of demineralized bone extract to induce endochondral osteogenesis in vivo in an extraskelatal site. This particular family member plays an important role in the onset of endochondral bone formation in humans, and a reduction in expression has been associated with a variety of bone diseases, including the heritable disorder Fibrodysplasia Ossificans Progressiva. Alternative splicing in the 5' untranslated region of this gene has been described and three variants are described, all encoding an identical protein.


研究领域 (Research Area)

TGF-beta signaling pathway

图片 (Image Data)




Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);




Western blot analysis using BMP4 mouse mAb against A549 (1) and C6 (2) cell lysate.

产品名称: BMP4 Mouse Monoclonal Antibody
产品货号: AMM81584



Immunofluorescence analysis of HeLa cells using BMP4 mouse mAb. Blue: DRAQ5 fluorescent DNA dye. Red: Actin filaments have been labeled with Alexa Fluor- 555 phalloidin.



Immunofluorescence analysis of HeLa cells using BMP4 mouse mAb (green). Red: Actin filaments have been labeled with Alexa Fluor- 555 phalloidin. Secondary antibody from Fisher (Cat#: 35503)

注意事项 (Note)

For research use only .