

产品名称: FOXM1 Mouse Monoclonal Antibody
产品货号: AMM81593



产品概述 (Summary)

产品名称 (Production Name)	FOXM1 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	ELISA,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	FOXM1
别名 (Alternative Names)	MPP2; TGT3; HFH11; HNF-3; INS-1; MPP-2; PIG29; FKHL16; FOXM1B; HFH-11; TRIDENT; MPHOSPH2
基因 ID (Gene ID)	2305.0
蛋白 ID (SwissProt ID)	Q08050.Purified recombinant fragment of human FOXM1 (AA: 649-748) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	84.3kDa

研究背景 (Background)


产品名称: FOXM1 Mouse Monoclonal Antibody
产品货号: AMM81593




The protein encoded by this gene is a transcriptional activator involved in cell proliferation. The encoded protein is phosphorylated in M phase and regulates the expression of several cell cycle genes, such as cyclin B1 and cyclin D1. Several transcript variants encoding different isoforms have been found for this gene.

研究领域 (Research Area)

图片 (Image Data)



Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);



Flow cytometric analysis of HeLa cells using FOX M1 mouse mAb (green) and negative control (red).

注意事项 (Note)

For research use only .