

产品名称: ATG14L Mouse Monoclonal Antibody
产品货号: AMM81601



产品概述 (Summary)

产品名称 (Production Name)	ATG14L Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	IHC,ELISA
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ATG14L
别名 (Alternative Names)	ATG14; BARKOR; KIAA0831
基因 ID (Gene ID)	22863.0
蛋白 ID (SwissProt ID)	Q6ZNE5.Purified recombinant fragment of human ATG14L (AA: 43-303) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:200-1:1000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	55.3kDa

研究背景 (Background)

ATG14 (Autophagy Related 14) is a Protein Coding gene. Diseases associated with ATG14 include human granulocytic

产品名称: ATG14L Mouse Monoclonal Antibody
产品货号: AMM81601




anaplasmosis. Among its related pathways are Senescence and Autophagy and Regulation of autophagy.


研究领域 (Research Area)

Autophagy

图片 (Image Data)



Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);



Immunohistochemical analysis of paraffin-embedded human testis tissues using ATG14L mouse mAb with DAB staining.

注意事项 (Note)

For research use only .