

产品名称: IL28A Mouse Monoclonal Antibody
产品货号: AMM81609



产品概述 (Summary)

产品名称 (Production Name)	IL28A Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IHC,ELISA
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	IL28A
别名 (Alternative Names)	IFNL2; IL-28A
基因 ID (Gene ID)	282616.0
蛋白 ID (SwissProt ID)	Q8IZJ0. Purified recombinant fragment of human IL28A (AA: 1-200) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:200-1:1000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	22.3kDa

研究背景 (Background)

This gene encodes a cytokine distantly related to type I interferons and the IL-10 family. This gene, interleukin 28B (IL28B),

产品名称: IL28A Mouse Monoclonal Antibody
产品货号: AMM81609




and interleukin 29 (IL29) are three closely related cytokine genes that form a cytokine gene cluster on a chromosomal region mapped to 19q13. Expression of the cytokines encoded by the three genes can be induced by viral infection. All three cytokines have been shown to interact with a heterodimeric class II cytokine receptor that consists of interleukin 10 receptor, beta (IL10RB) and interleukin 28 receptor, alpha (IL28RA).


研究领域 (Research Area)

Jak-STAT signaling pathway


图片 (Image Data)



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)




Western blot analysis using IL28A mouse mAb against SK-Br-3 (1) cell lysate.



产品名称: IL28A Mouse Monoclonal Antibody
产品货号: AMM81609



Immunohistochemical analysis of paraffin-embedded human brain tissues using IL28A mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human esophageal cancer tissues using IL28A mouse mAb with DAB staining.

注意事项 (Note)

For research use only .