

产品名称: KPNA2 Mouse Monoclonal Antibody
产品货号: AMM81616



产品概述 (Summary)

产品名称 (Production Name)	KPNA2 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IHC,ICC,ELISA,FC
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Phosphorylated
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	KPNA2
别名 (Alternative Names)	QIP2; RCH1; IPOA1; SRP1alpha
基因 ID (Gene ID)	3838.0
蛋白 ID (SwissProt ID)	P52292 .Purified recombinant fragment of human KPNA2 (AA: 1-530) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:200-1:1000,ICC 1:200-1:1000,ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	58kDa

研究背景 (Background)


产品名称: KPNA2 Mouse Monoclonal Antibody
产品货号: AMM81616




The import of proteins into the nucleus is a process that involves at least 2 steps. The first is an energy-independent docking of the protein to the nuclear envelope and the second is an energy-dependent translocation through the nuclear pore complex. Imported proteins require a nuclear localization sequence (NLS) which generally consists of a short region of basic amino acids or 2 such regions spaced about 10 amino acids apart. Proteins involved in the first step of nuclear import have been identified in different systems. These include the *Xenopus* protein importin and its yeast homolog, SRP1 (a suppressor of certain temperature-sensitive mutations of RNA polymerase I in *Saccharomyces cerevisiae*), which bind to the NLS. KPNA2 protein interacts with the NLSs of DNA helicase Q1 and SV40 T antigen and may be involved in the nuclear transport of proteins. KPNA2 also may play a role in V(D)J recombination.

研究领域 (Research Area)

图片 (Image Data)




Black line: Control Antigen (100 ng); Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng)




Western blot analysis using KPNA2 mouse mAb against Hela (1), HEK293 (2), and NIH/3T3 (3) cell lysate.


产品名称: KPNA2 Mouse Monoclonal Antibody
产品货号: AMM81616



Immunofluorescence analysis of Hela cells using KPNA2 mouse mAb (green). Blue: DRAQ5 fluorescent DNA dye. Red: Actin filaments have been labeled with Alexa Fluor- 555 phalloidin. Secondary antibody from Fisher (Cat#: 35503)




Flow cytometric analysis of Hela cells using KPNA2 mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human lung cancer tissues using KPNA2 mouse mAb with DAB staining.

产品名称: KPNA2 Mouse Monoclonal Antibody
产品货号: AMM81616



Immunohistochemical analysis of paraffin-embedded human esophageal cancer tissues using KPNA2 mouse mAb with DAB staining.

注意事项 (Note)

For research use only .