

产品名称: KDM6A Mouse Monoclonal Antibody
产品货号: AMM81624



产品概述 (Summary)

产品名称 (Production Name)	KDM6A Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG2a
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	KDM6A
别名 (Alternative Names)	UTX; KABUK2; bA386N14.2
基因 ID (Gene ID)	7403.0
蛋白 ID (SwissProt ID)	O15550.Purified recombinant fragment of human KDM6A (AA: 1252-1401) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	154kDa

研究背景 (Background)

This gene is located on the X chromosome and is the corresponding locus to a Y-linked gene which encodes a


产品名称: KDM6A Mouse Monoclonal Antibody
产品货号: AMM81624




tetratricopeptide repeat (TPR) protein. The encoded protein of this gene contains a JmjC-domain and catalyzes the demethylation of tri/dimethylated histone H3. Multiple alternatively spliced transcript variants have been found for this gene.

研究领域 (Research Area)

图片 (Image Data)



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Western blot analysis using KDM6A mouse mAb against HepG2 (1) and C6 (2) cell lysate.

注意事项 (Note)

For research use only .