

产品名称: OTUD4 Mouse Monoclonal Antibody
产品货号: AMM81634



产品概述 (Summary)

产品名称 (Production Name)	OTUD4 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,ELISA,FC
种属反应性 (Reactivity)	Human,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	OTUD4
别名 (Alternative Names)	HIN1; DUBA6; HSHIN1
基因 ID (Gene ID)	54726.0
蛋白 ID (SwissProt ID)	Q01804.Purified recombinant fragment of human OTUD4 (AA: 815-1049) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	124kDa

研究背景 (Background)

Alternatively spliced transcript variants have been found for this gene. The smaller protein isoform encoded by the shorter


产品名称: OTUD4 Mouse Monoclonal Antibody
产品货号: AMM81634




transcript variant is found only in HIV-1 infected cells.

研究领域 (Research Area)


图片 (Image Data)



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Western blot analysis using OTUD4 mouse mAb against HepG2 (1) and C6 (2) cell lysate.



Flow cytometric analysis of Hela cells using OTUD4 mouse mAb (green) and negative control (red).

产品名称: OTUD4 Mouse Monoclonal Antibody
产品货号: AMM81634



注意事项 (Note)

For research use only .