

产品名称: Rab1b Mouse Monoclonal Antibody
产品货号: AMM81635



产品概述 (Summary)

产品名称 (Production Name)	Rab1b Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	IHC,ELISA,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG2b
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Rab1b
别名 (Alternative Names)	Rab1b
基因 ID (Gene ID)	81876.0
蛋白 ID (SwissProt ID)	Q9H0U4.Purified recombinant fragment of human Rab1b (AA: 60-201) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:200-1:1000,ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	22.2kDa

研究背景 (Background)

Members of the RAB protein family, such as RAB1B, are low molecular mass monomeric GTPases localized on the

产品名称: Rab1b Mouse Monoclonal Antibody
产品货号: AMM81635




cytoplasmic surfaces of distinct membrane-bound organelles. RAB1B functions in the early secretory pathway and is essential for vesicle transport between the endoplasmic reticulum (ER) and Golgi


研究领域 (Research Area)

Autophagy


图片 (Image Data)



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)




Flow cytometric analysis of Hela cells using Rab1b mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human bladder cancer tissues using Rab1b mouse mAb with DAB staining.

产品名称: Rab1b Mouse Monoclonal Antibody
产品货号: AMM81635



Immunohistochemical analysis of paraffin-embedded human esophageal cancer tissues using Rab1b mouse mAb with DAB staining.

注意事项 (Note)

For research use only .