

产品名称: DDX5 Mouse Monoclonal Antibody
产品货号: AMM81651



产品概述 (Summary)

产品名称 (Production Name)	DDX5 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IHC,ELISA,FC
种属反应性 (Reactivity)	Human,Mouse,Monkey

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG2a
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	DDX5
别名 (Alternative Names)	p68; HLR1; G17P1; HUMP68
基因 ID (Gene ID)	1655.0
蛋白 ID (SwissProt ID)	P17844.Purified recombinant fragment of human DDX5 (AA: 475-614) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:200-1:1000,ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	69.1kDa

研究背景 (Background)

DEAD box proteins, characterized by the conserved motif Asp-Glu-Ala-Asp (DEAD), are putative RNA helicases. They are


产品名称: DDX5 Mouse Monoclonal Antibody
产品货号: AMM81651




implicated in a number of cellular processes involving alteration of RNA secondary structure, such as translation initiation, nuclear and mitochondrial splicing, and ribosome and spliceosome assembly. Based on their distribution patterns, some members of this family are believed to be involved in embryogenesis, spermatogenesis, and cellular growth and division. This gene encodes a DEAD box protein, which is a RNA-dependent ATPase, and also a proliferation-associated nuclear antigen, specifically reacting with the simian virus 40 tumor antigen. This gene consists of 13 exons, and alternatively spliced transcripts containing several intron sequences have been detected, but no isoforms encoded by these transcripts have been identified.

研究领域 (Research Area)

图片 (Image Data)




Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)




Western blot analysis using DDX5 mouse mAb against HT-29 (1), Hela (2), NIH/3T3 (3), COS7 (4), SW620 (5), Jurkat (6), A431 (7), and MCF-7 (8) cell lysate.

产品名称: DDX5 Mouse Monoclonal Antibody
产品货号: AMM81651



Flow cytometric analysis of Hela cells using DDX5 mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human bladder cancer tissues using DDX5 mouse mAb with DAB staining.

注意事项 (Note)

For research use only .