

**产品名称: PELP1 Mouse Monoclonal Antibody**  
**产品货号: AMM81662**



## 产品概述 (Summary)

|                        |                                 |
|------------------------|---------------------------------|
| 产品名称 (Production Name) | PELP1 Mouse Monoclonal Antibody |
| 描述 (Description)       | Mouse monoclonal Antibody       |
| 宿主 (Host)              | Mouse                           |
| 应用 (Application)       | ELISA,FC                        |
| 种属反应性 (Reactivity)     | Human                           |

## 产品性能 (Performance)

|                     |                                                                                          |
|---------------------|------------------------------------------------------------------------------------------|
| 偶联物 (Conjugation)   | Unconjugated                                                                             |
| 修饰 (Modification)   | Unmodified                                                                               |
| 同种型 (Isotype)       | Mouse IgG2a                                                                              |
| 克隆 (Clonality)      | Monoclonal                                                                               |
| 形式 (Form)           | Liquid                                                                                   |
| 存放说明 (Storage)      | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles. |
| 储存溶液 (Buffer)       | Purified antibody in PBS with 0.05% sodium azide                                         |
| 纯化方式 (Purification) | Affinity Purification                                                                    |

## 免疫原信息 (Immunogen)

|                        |                                                                                           |
|------------------------|-------------------------------------------------------------------------------------------|
| 基因名 (Gene Name)        | PELP1                                                                                     |
| 别名 (Alternative Names) | MNAR; P160                                                                                |
| 基因 ID (Gene ID)        | 27043.0                                                                                   |
| 蛋白 ID (SwissProt ID)   | Q8IZL8.Purified recombinant fragment of human PELP1 (AA: 1031-1180) expressed in E. Coli. |

## 产品应用 (Application)

|                          |                                     |
|--------------------------|-------------------------------------|
| 稀释比 (Dilution Ratio)     | ELISA 1:5000-1:20000,FC 1:200-1:400 |
| 蛋白分子量 (Molecular Weight) | 120kDa                              |

## 研究背景 (Background)

This gene encodes a transcription factor which coactivates transcription of estrogen receptor responsive genes and


产品名称: PELP1 Mouse Monoclonal Antibody  
产品货号: AMM81662




corepresses genes activated by other hormone receptors or sequence-specific transcription factors. Expression of this gene is regulated by both members of the estrogen receptor family. This gene may be involved in the progression of several types of cancer. Alternative splicing results in multiple transcript variants.

## 研究领域 (Research Area)

## 图片 (Image Data)



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Flow cytometric analysis of Hela cells using PELP1 mouse mAb (green) and negative control (red).

## 注意事项 (Note)

For research use only .