

产品名称: PELP1 Mouse Monoclonal Antibody
产品货号: AMM81663



产品概述 (Summary)

产品名称 (Production Name)	PELP1 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	ELISA,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG2b
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PELP1
别名 (Alternative Names)	MNAR; P160
基因 ID (Gene ID)	27043.0
蛋白 ID (SwissProt ID)	Q8IZL8.Purified recombinant fragment of human PELP1 (AA: 1031-1180) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	120kDa

研究背景 (Background)

This gene encodes a transcription factor which coactivates transcription of estrogen receptor responsive genes and


产品名称: PELP1 Mouse Monoclonal Antibody
产品货号: AMM81663




corepresses genes activated by other hormone receptors or sequence-specific transcription factors. Expression of this gene is regulated by both members of the estrogen receptor family. This gene may be involved in the progression of several types of cancer. Alternative splicing results in multiple transcript variants.

研究领域 (Research Area)

图片 (Image Data)



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Flow cytometric analysis of Hela cells using PELP1 mouse mAb (green) and negative control (red).

注意事项 (Note)

For research use only .