

产品名称: XRN2 Mouse Monoclonal Antibody
产品货号: AMM81673



产品概述 (Summary)

产品名称 (Production Name)	XRN2 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IHC,ELISA
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	XRN2
别名 (Alternative Names)	XRN2
基因 ID (Gene ID)	22803.0
蛋白 ID (SwissProt ID)	Q9H0D6.Purified recombinant fragment of human XRN2 (AA: 398-547) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:200-1:1000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	108.5kDa

研究背景 (Background)

This gene encodes a 5'-3' exonuclease that promotes transcription termination at cotranscriptional cleavage sites.


产品名称: XRN2 Mouse Monoclonal Antibody
产品货号: AMM81673




Alternative splicing results in multiple transcript variants encoding different isoforms.

研究领域 (Research Area)


图片 (Image Data)



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)




Western blot analysis using XRN2 mouse mAb against Raw264.7 (1), HEK293 (2), NTERA-2 (3), LNCap (4), HepG2 (5), HEK293 (6), and Hela (7) cell lysate.



Immunohistochemical analysis of paraffin-embedded human colon cancer tissues using XRN2 mouse mAb with DAB staining.

产品名称: XRN2 Mouse Monoclonal Antibody
产品货号: AMM81673



Immunohistochemical analysis of paraffin-embedded human ovarian cancer tissues using XRN2 mouse mAb with DAB staining.

注意事项 (Note)

For research use only .