

产品名称: ARF1 Mouse Monoclonal Antibody
产品货号: AMM81679



产品概述 (Summary)

产品名称 (Production Name)	ARF1 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	IHC, ICC, ELISA, FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG2a
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ARF1
别名 (Alternative Names)	ARF1
基因 ID (Gene ID)	375.0
蛋白 ID (SwissProt ID)	P84077. Purified recombinant fragment of human ARF1 (AA: 76-182) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:200-1:1000, ICC 1:200-1:1000, ELISA 1:5000-1:20000, FC 1:200-1:400
蛋白分子量 (Molecular Weight)	20.7kDa

研究背景 (Background)

ADP-ribosylation factor 1 (ARF1) is a member of the human ARF gene family. The family members encode small guanine


产品名称: ARF1 Mouse Monoclonal Antibody
产品货号: AMM81679




nucleotide-binding proteins that stimulate the ADP-ribosyltransferase activity of cholera toxin and play a role in vesicular trafficking as activators of phospholipase D. The gene products, including 6 ARF proteins and 11 ARF-like proteins, constitute a family of the RAS superfamily. The ARF proteins are categorized as class I (ARF1, ARF2 and ARF3), class II (ARF4 and ARF5) and class III (ARF6), and members of each class share a common gene organization. The ARF1 protein is localized to the Golgi apparatus and has a central role in intra-Golgi transport. Multiple alternatively spliced transcript variants encoding the same protein have been found for this gene.

研究领域 (Research Area)

图片 (Image Data)




Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)




Immunofluorescence analysis of HL-7702 cells using ARF1 mouse mAb (green). Blue: DRAQ5 fluorescent DNA dye. Red: Actin filaments have been labeled with Alexa Fluor- 555 phalloidin. Secondary antibody from Fisher (Cat#: 35503)


产品名称: ARF1 Mouse Monoclonal Antibody
产品货号: AMM81679



Immunofluorescence analysis of SK-OV-3 cells using ARF1 mouse mAb (green). Blue: DRAQ5 fluorescent DNA dye. Red: Actin filaments have been labeled with Alexa Fluor- 555 phalloidin. Secondary antibody from Fisher (Cat#: 35503)




Flow cytometric analysis of Hela cells using ARF1 mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human bladder cancer tissues using ARF1 mouse mAb with DAB staining.

产品名称: ARF1 Mouse Monoclonal Antibody
产品货号: AMM81679



Immunohistochemical analysis of paraffin-embedded human stomach cancer tissues using ARF1 mouse mAb with DAB staining.

注意事项 (Note)

For research use only .