

产品名称: B3GAT1 Mouse Monoclonal Antibody
产品货号: AMM81713



产品概述 (Summary)

产品名称 (Production Name)	B3GAT1 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,ELISA,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	B3GAT1
别名 (Alternative Names)	NK1; CD57; HNK1; LEU7; NK-1; GLCATP; GLCUATP
基因 ID (Gene ID)	27087.0
蛋白 ID (SwissProt ID)	Q9P2W7.Purified recombinant fragment of human B3GAT1 (AA: 193-334) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	38.2kDa

研究背景 (Background)

The protein encoded by this gene is a member of the glucuronyltransferase gene family. These enzymes exhibit strict


产品名称: B3GAT1 Mouse Monoclonal Antibody
产品货号: AMM81713




acceptor specificity, recognizing nonreducing terminal sugars and their anomeric linkages. This gene product functions as the key enzyme in a glucuronyl transfer reaction during the biosynthesis of the carbohydrate epitope HNK-1 (human natural killer-1, also known as CD57 and LEU7). Alternate transcriptional splice variants have been characterized.

研究领域 (Research Area)


图片 (Image Data)



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Western blot analysis using B3GAT1 mouse mAb against MOLT4 (1), Raji (2), and Ramos (3) cell lysate.



Flow cytometric analysis of Hela cells using B3GAT1 mouse mAb (green) and negative control (red).

产品名称: B3GAT1 Mouse Monoclonal Antibody
产品货号: AMM81713



注意事项 (Note)

For research use only .