

产品名称: NR1I2 Mouse Monoclonal Antibody
产品货号: AMM81773



产品概述 (Summary)

产品名称 (Production Name)	NR1I2 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	ICC,ELISA,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG2b
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	NR1I2
别名 (Alternative Names)	BXR; PAR; PRR; PXR; SAR; SXR; ONR1; PAR1; PAR2; PARq
基因 ID (Gene ID)	8856.0
蛋白 ID (SwissProt ID)	O75469.Purified recombinant fragment of human NR1I2 (AA: 1-142) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	ICC 1:200-1:1000,ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	49.8kDa

研究背景 (Background)

This gene product belongs to the nuclear receptor superfamily, members of which are transcription factors characterized by


产品名称: NR1I2 Mouse Monoclonal Antibody
产品货号: AMM81773




a ligand-binding domain and a DNA-binding domain. The encoded protein is a transcriptional regulator of the cytochrome P450 gene CYP3A4, binding to the response element of the CYP3A4 promoter as a heterodimer with the 9-cis retinoic acid receptor RXR. It is activated by a range of compounds that induce CYP3A4, including dexamethasone and rifampicin. Several alternatively spliced transcripts encoding different isoforms, some of which use non-AUG (CUG) translation initiation codon, have been described for this gene. Additional transcript variants exist, however, they have not been fully characterized.

研究领域 (Research Area)

图片 (Image Data)




Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Immunofluorescence analysis of HeLa cells using NR1I2 mouse mAb (green). Blue: DRAQ5 fluorescent DNA dye. Red: Actin filaments have been labeled with Alexa Fluor- 555 phalloidin. Secondary antibody from Fisher (Cat#: 35503)

产品名称: NR1I2 Mouse Monoclonal Antibody
产品货号: AMM81773



Flow cytometric analysis of HepG2 cells using NR1I2 mouse mAb (green) and negative control (red).

注意事项 (Note)

For research use only .