

产品名称: CPV-VP2 Mouse Monoclonal Antibody
产品货号: AMM81775



产品概述 (Summary)

产品名称 (Production Name)	CPV-VP2 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IHC,ICC,ELISA,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CPV-VP2
别名 (Alternative Names)	Canine Parvovirus VP2
基因 ID (Gene ID)	EF011664.1
蛋白 ID (SwissProt ID)	.Purified recombinant fragment of Canine Parvovirus CPV-VP2 (AA: 1-128) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:200-1:1000,ICC 1:200-1:1000,ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	64.6kDa


研究背景 (Background)

产品名称: CPV-VP2 Mouse Monoclonal Antibody
产品货号: AMM81775




研究领域 (Research Area)


图片 (Image Data)



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)




Western blot analysis using CPV-VP2 mouse mAb against SW480 (1), A431 (2), K562 (3), and BCBL-1 (4) cell lysate.




Immunofluorescence analysis of HeLa cells using CPV-VP2 mouse mAb (green). Blue: DRAQ5 fluorescent DNA dye. Red: Actin filaments have been labeled with Alexa Fluor- 555 phalloidin. Secondary antibody from Fisher (Cat#: 35503)


产品名称: CPV-VP2 Mouse Monoclonal Antibody
产品货号: AMM81775



Flow cytometric analysis of Hela cells using CPV-VP2 mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human cervical cancer tissues using CPV-VP2 mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human colon cancer tissues using CPV-VP2 mouse mAb with DAB staining.

注意事项 (Note)

For research use only .