

产品名称: RALB Mouse Monoclonal Antibody
产品货号: AMM81777



产品概述 (Summary)

产品名称 (Production Name)	RALB Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IHC,ELISA,FC
种属反应性 (Reactivity)	Human,Mouse,Monkey

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG2b
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	RALB
别名 (Alternative Names)	RALB
基因 ID (Gene ID)	5899.0
蛋白 ID (SwissProt ID)	P11234.Purified recombinant fragment of human RALB (AA: 89-206) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:200-1:1000,ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	23.4kDa

研究背景 (Background)

This gene encodes a GTP-binding protein that belongs to the small GTPase superfamily and Ras family of proteins. GTP-

产品名称: RALB Mouse Monoclonal Antibody
产品货号: AMM81777




binding proteins mediate the transmembrane signaling initiated by the occupancy of certain cell surface receptors.


研究领域 (Research Area)

Apoptosis


图片 (Image Data)



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)




Western blot analysis using RALB mouse mAb against HepG2 (1), COS7 (2), NIH/3T3 (3), A549 (4), U251 (5), HT-29 (6), HEK293 (7), and LOVO (8) cell lysate.




Flow cytometric analysis of Hela cells using RALB mouse mAb (green) and negative control (red).

产品名称: RALB Mouse Monoclonal Antibody
产品货号: AMM81777



Immunohistochemical analysis of paraffin-embedded human cervical cancer tissues using RALB mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human bladder cancer tissues using RALB mouse mAb with DAB staining.

注意事项 (Note)

For research use only .