

产品名称: AFAP1L2 Mouse Monoclonal Antibody
产品货号: AMM81778



产品概述 (Summary)

产品名称 (Production Name)	AFAP1L2 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,ELISA,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	AFAP1L2
别名 (Alternative Names)	XB130; KIAA1914; CTB-1144G6.4
基因 ID (Gene ID)	84632.0
蛋白 ID (SwissProt ID)	Q8N4X5 .Purified recombinant fragment of human AFAP1L2 (AA: 674-818) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	91.3kDa

研究背景 (Background)

AFAP1L2 (Actin Filament Associated Protein 1 Like 2) is a Protein Coding gene. GO annotations related to this gene include


产品名称: AFAP1L2 Mouse Monoclonal Antibody
产品货号: AMM81778




SH3 domain binding and protein tyrosine kinase activator activity. An important paralog of this gene is AFAP1.

研究领域 (Research Area)


图片 (Image Data)



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Western blot analysis using AFAP1L2 mouse mAb against A549 (1) cell lysate.



Flow cytometric analysis of Hela cells using AFAP1L2 mouse mAb (green) and negative control (red).

产品名称: AFAP1L2 Mouse Monoclonal Antibody
产品货号: AMM81778



注意事项 (Note)

For research use only .