

产品名称: ADRB2 Mouse Monoclonal Antibody
产品货号: AMM81790



产品概述 (Summary)

产品名称 (Production Name)	ADRB2 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,ELISA,FC
种属反应性 (Reactivity)	Human,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ADRB2
别名 (Alternative Names)	BAR; B2AR; ADRBR; ADRB2R; BETA2AR
基因 ID (Gene ID)	154.0
蛋白 ID (SwissProt ID)	P07550.Purified recombinant fragment of human ADRB2 (AA: 302-413) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	46.5kDa

研究背景 (Background)

This gene encodes beta-2-adrenergic receptor which is a member of the G protein-coupled receptor superfamily. This


产品名称: ADRB2 Mouse Monoclonal Antibody
产品货号: AMM81790




receptor is directly associated with one of its ultimate effectors, the class C L-type calcium channel Ca(V)1.2. This receptor-channel complex also contains a G protein, an adenylyl cyclase, cAMP-dependent kinase, and the counterbalancing phosphatase, PP2A. The assembly of the signaling complex provides a mechanism that ensures specific and rapid signaling by this G protein-coupled receptor. This gene is intronless. Different polymorphic forms, point mutations, and/or downregulation of this gene are associated with nocturnal asthma, obesity and type 2 diabetes.

研究领域 (Research Area)

图片 (Image Data)




Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Western blot analysis using ADRB2 mouse mAb against C6 (1) cell lysate.

产品名称: ADRB2 Mouse Monoclonal Antibody
产品货号: AMM81790



Flow cytometric analysis of MCF-7 cells using ADRB2 mouse mAb (green) and negative control (red).

注意事项 (Note)

For research use only .