

产品名称: NEFM Mouse Monoclonal Antibody
产品货号: AMM81800



产品概述 (Summary)

产品名称 (Production Name)	NEFM Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	ICC,ELISA,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	NEFM
别名 (Alternative Names)	NFM; NEF3; NF-M
基因 ID (Gene ID)	4741.0
蛋白 ID (SwissProt ID)	P07197.Purified recombinant fragment of human NEFM (AA: 779-916) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	ICC 1:50-1:200,ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	102.5kDa

研究背景 (Background)

Neurofilaments are type IV intermediate filament heteropolymers composed of light, medium, and heavy chains.


产品名称: NEFM Mouse Monoclonal Antibody
产品货号: AMM81800




Neurofilaments comprise the axoskeleton and functionally maintain neuronal caliber. They may also play a role in intracellular transport to axons and dendrites. This gene encodes the medium neurofilament protein. This protein is commonly used as a biomarker of neuronal damage. Alternative splicing results in multiple transcript variants encoding distinct isoforms.

研究领域 (Research Area)

图片 (Image Data)




Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)




Immunofluorescence analysis of HeLa cells using NEFM mouse mAb (green). Blue: DRAQ5 fluorescent DNA dye. Red: Actin filaments have been labeled with Alexa Fluor- 555 phalloidin. Secondary antibody from Fisher (Cat#: 35503)

产品名称: NEFM Mouse Monoclonal Antibody
产品货号: AMM81800



Flow cytometric analysis of Hela cells using NEFM mouse mAb (green) and negative control (red).



Flow cytometric analysis of Raji cells using NEFM mouse mAb (green) and negative control (red).

注意事项 (Note)

For research use only .