

产品名称: SAG2 Mouse Monoclonal Antibody
产品货号: AMM81801



产品概述 (Summary)

产品名称 (Production Name)	SAG2 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,ELISA,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	SAG2
别名 (Alternative Names)	Toxoplasma gondii surface antigen 2B
基因 ID (Gene ID)	KC607825.1
蛋白 ID (SwissProt ID)	.Purified recombinant fragment of Toxoplasma gondii surface antigen 2B SAG2 (AA: 1-130) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	13.4kDa


研究背景 (Background)

产品名称: SAG2 Mouse Monoclonal Antibody
产品货号: AMM81801




研究领域 (Research Area)


图片 (Image Data)



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Western blot analysis using SAG2 mouse mAb against SW480 (1), A431 (2), and K562 (3) cell lysate.



Flow cytometric analysis of HeLa cells using SAG2 mouse mAb (green) and negative control (red).

注意事项 (Note)

产品名称: SAG2 Mouse Monoclonal Antibody
产品货号: AMM81801



For research use only .