

产品名称: ASH2L Mouse Monoclonal Antibody
产品货号: AMM81819



产品概述 (Summary)

产品名称 (Production Name)	ASH2L Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IHC,ELISA,FC
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ASH2L
别名 (Alternative Names)	ASH2; Bre2; ASH2L1; ASH2L2
基因 ID (Gene ID)	9070.0
蛋白 ID (SwissProt ID)	Q9UBL3.Purified recombinant fragment of human ASH2L (AA: 493-628) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:200-1:1000,ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	68.7kDa

研究背景 (Background)

ASH2L (Ash2 (Absent, Small, Or Homeotic)-Like (Drosophila)) is a Protein Coding gene. Diseases associated with ASH2L


产品名称: ASH2L Mouse Monoclonal Antibody
产品货号: AMM81819




include Kabuki Syndrome 1. Among its related pathways are Developmental Biology and Signaling by Wnt. GO annotations related to this gene include transcription regulatory region DNA binding and histone methyltransferase activity (H3-K4 specific).

研究领域 (Research Area)


图片 (Image Data)



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)




Western blot analysis using ASH2L mouse mAb against K562 (1) and F9 (2) cell lysate.




Flow cytometric analysis of K562 cells using ASH2L mouse mAb (green) and negative control (red).

产品名称: ASH2L Mouse Monoclonal Antibody
产品货号: AMM81819



Immunohistochemical analysis of paraffin-embedded human colon cancer tissues using ASH2L mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human rectum cancer tissues using ASH2L mouse mAb with DAB staining.

注意事项 (Note)

For research use only .