

产品名称: ASF1B Mouse Monoclonal Antibody
产品货号: AMM81843



产品概述 (Summary)

产品名称 (Production Name)	ASF1B Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,ELISA,FC
种属反应性 (Reactivity)	Human,Mouse,Monkey

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ASF1B
别名 (Alternative Names)	CIA-II
基因 ID (Gene ID)	55723.0
蛋白 ID (SwissProt ID)	Q9NVP2.Purified recombinant fragment of human ASF1B (AA: 1-202) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	22.4kDa

研究背景 (Background)

This gene encodes a member of the H3/H4 family of histone chaperone proteins and is similar to the anti-silencing


产品名称: ASF1B Mouse Monoclonal Antibody
产品货号: AMM81843




function-1 gene in yeast. The encoded protein is the substrate of the tousel-like kinase family of cell cycle-regulated kinases, and may play a key role in modulating the nucleosome structure of chromatin by ensuring a constant supply of histones at sites of nucleosome assembly.

研究领域 (Research Area)


图片 (Image Data)



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Western blot analysis using ASF1B mouse mAb against Hela (1), COS7 (2), HCT116 (3), and CHO3D10 (4) cell lysate.



Flow cytometric analysis of K562 cells using ASF1B mouse mAb (green) and negative control (red).

产品名称: ASF1B Mouse Monoclonal Antibody
产品货号: AMM81843



注意事项 (Note)

For research use only .