

产品名称: IGLC2 Mouse Monoclonal Antibody
产品货号: AMM81844



产品概述 (Summary)

产品名称 (Production Name)	IGLC2 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	IHC,ELISA
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	IGLC2
别名 (Alternative Names)	IGLC
基因 ID (Gene ID)	3538.0
蛋白 ID (SwissProt ID)	P0DOY2.Purified recombinant fragment of human IGLC2 (AA: 1-106) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:200-1:1000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	11.3kDa

研究背景 (Background)

IGLC2 (Immunoglobulin Lambda Constant 2) is a Protein Coding gene. Among its related pathways are Interleukin-3, 5 and


产品名称: IGLC2 Mouse Monoclonal Antibody
产品货号: AMM81844




GM-CSF signaling and Vesicle-mediated transport. GO annotations related to this gene include antigen binding. An important paralog of this gene is IGLC3.

研究领域 (Research Area)


图片 (Image Data)



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Immunohistochemical analysis of paraffin-embedded human colon tissues using IGLC2 mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human endometrial cancer tissues using IGLC2 mouse mAb with DAB staining.

产品名称: IGLC2 Mouse Monoclonal Antibody
产品货号: AMM81844



注意事项 (Note)

For research use only .