

产品名称: TFAP2C Mouse Monoclonal Antibody
产品货号: AMM81845



产品概述 (Summary)

产品名称 (Production Name)	TFAP2C Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IHC,ELISA
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TFAP2C
别名 (Alternative Names)	ERF1; TFAP2G; hAP-2g; AP2-GAMMA
基因 ID (Gene ID)	7022.0
蛋白 ID (SwissProt ID)	Q92754.Purified recombinant fragment of human TFAP2C (AA: 341-450) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:200-1:1000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	49.2kDa

研究背景 (Background)

The protein encoded by this gene is a sequence-specific DNA-binding transcription factor involved in the activation of


产品名称: TFAP2C Mouse Monoclonal Antibody
产品货号: AMM81845




several developmental genes. The encoded protein can act as either a homodimer or heterodimer with other family members and is induced during retinoic acid-mediated differentiation. It plays a role in the development of the eyes, face, body wall, limbs, and neural tube.

研究领域 (Research Area)


图片 (Image Data)



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Western blot analysis using TFAP2C mouse mAb against SK-Br-3 (1) cell lysate.



Immunohistochemical analysis of paraffin-embedded human esophageal cancer tissues using TFAP2C mouse mAb with DAB

产品名称: TFAP2C Mouse Monoclonal Antibody
产品货号: AMM81845



staining.

注意事项 (Note)

For research use only .