

产品名称: KDM4B Mouse Monoclonal Antibody
产品货号: AMM81847



产品概述 (Summary)

产品名称 (Production Name)	KDM4B Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	ELISA,FC
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	KDM4B
别名 (Alternative Names)	JMJD2B; TDRD14B
基因 ID (Gene ID)	23030.0
蛋白 ID (SwissProt ID)	O94953.Purified recombinant fragment of human KDM4B (AA: 960-1096) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	122kDa

研究背景 (Background)

KDM4B (Lysine Demethylase 4B) is a Protein Coding gene. Among its related pathways are Chromatin organization and


产品名称: KDM4B Mouse Monoclonal Antibody
产品货号: AMM81847




DNA Double-Strand Break Repair. GO annotations related to this gene include dioxygenase activity. An important paralog of this gene is KDM4C.

研究领域 (Research Area)

图片 (Image Data)



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Flow cytometric analysis of NIH/3T3 cells using KDM4B mouse mAb (green) and negative control (red).

注意事项 (Note)

For research use only .