

产品名称: RIPK4 Mouse Monoclonal Antibody
产品货号: AMM81853



产品概述 (Summary)

产品名称 (Production Name)	RIPK4 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	ELISA,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	RIPK4
别名 (Alternative Names)	DIK; PKK; PPS2; RIP4; ANKK2; NKRD3; ANKRD3
基因 ID (Gene ID)	54101.0
蛋白 ID (SwissProt ID)	Q9H4D1.Purified recombinant fragment of human RIPK4 (AA: 675-784) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	91.6kDa

研究背景 (Background)

The protein encoded by this gene is a serine/threonine protein kinase that interacts with protein kinase C-delta. The


产品名称: RIPK4 Mouse Monoclonal Antibody
产品货号: AMM81853




encoded protein can also activate NFkappaB and is required for keratinocyte differentiation. This kinase undergoes autophosphorylation.

研究领域 (Research Area)

图片 (Image Data)



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Flow cytometric analysis of SK-N-SH cells using RIPK4 mouse mAb (green) and negative control (red).

注意事项 (Note)

For research use only .