

产品名称: ACVR1 Mouse Monoclonal Antibody
产品货号: AMM81888



产品概述 (Summary)

产品名称 (Production Name)	ACVR1 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	IHC, ICC, ELISA, FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ACVR1
别名 (Alternative Names)	FOP; ALK2; SKR1; TSRI; ACTRI; ACVR1A; ACVRLK2
基因 ID (Gene ID)	90.0
蛋白 ID (SwissProt ID)	Q04771. Purified recombinant fragment of human ACVR1 (AA: 21-120) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:200-1:1000, ICC 1:200-1:1000, ELISA 1:5000-1:20000, FC 1:200-1:400
蛋白分子量 (Molecular Weight)	57.2kDa

研究背景 (Background)

Activins are dimeric growth and differentiation factors which belong to the transforming growth factor-beta (TGF-beta)


产品名称: ACVR1 Mouse Monoclonal Antibody
产品货号: AMM81888




superfamily of structurally related signaling proteins. Activins signal through a heteromeric complex of receptor serine kinases which include at least two type I (I and IB) and two type II (II and IIB) receptors. These receptors are all transmembrane proteins, composed of a ligand-binding extracellular domain with cysteine-rich region, a transmembrane domain, and a cytoplasmic domain with predicted serine/threonine specificity. Type I receptors are essential for signaling; and type II receptors are required for binding ligands and for expression of type I receptors. Type I and II receptors form a stable complex after ligand binding, resulting in phosphorylation of type I receptors by type II receptors. This gene encodes activin A type I receptor which signals a particular transcriptional response in concert with activin type II receptors. Mutations in this gene are associated with fibrodysplasia ossificans progressive.

研究领域 (Research Area)

图片 (Image Data)




Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)




Immunofluorescence analysis of HeLa cells using ACVR1 mouse mAb. Blue: DRAQ5 fluorescent DNA dye. Red: Actin filaments have been labeled with Alexa Fluor- 555 phalloidin.


产品名称: ACVR1 Mouse Monoclonal Antibody
产品货号: AMM81888



Immunofluorescence analysis of HeLa cells using ACVR1 mouse mAb (green). Blue: DRAQ5 fluorescent DNA dye. Red: Actin filaments have been labeled with Alexa Fluor- 555 phalloidin. Secondary antibody from Fisher (Cat#: 35503)




Flow cytometric analysis of HeLa cells using ACVR1 mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human ovarian cancer tissues using ACVR1 mouse mAb with DAB staining.

产品名称: ACVR1 Mouse Monoclonal Antibody
产品货号: AMM81888



Immunohistochemical analysis of paraffin-embedded human endometrial cancer tissues using ACVR1 mouse mAb with DAB staining.

注意事项 (Note)

For research use only .