

产品名称: ATP1F1 Mouse Monoclonal Antibody
产品货号: AMM81910



产品概述 (Summary)

产品名称 (Production Name)	ATPIF1 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	IHC, ICC, ELISA, FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG2b
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ATPIF1
别名 (Alternative Names)	IP; ATP1; ATP1P
基因 ID (Gene ID)	93974.0
蛋白 ID (SwissProt ID)	Q9UII2. Purified recombinant fragment of human ATP1F1 (AA: 1-106) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:200-1:1000, ICC 1:100-1:500, ELISA 1:5000-1:20000, FC 1:200-1:400
蛋白分子量 (Molecular Weight)	12.3kDa

研究背景 (Background)

This gene encodes a mitochondrial ATPase inhibitor. Alternative splicing occurs at this locus and three transcript variants


产品名称: ATPIF1 Mouse Monoclonal Antibody
产品货号: AMM81910




encoding distinct isoforms have been identified.

研究领域 (Research Area)


图片 (Image Data)



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)




Immunofluorescence analysis of HeLa cells using ATPIF1 mouse mAb (green). Blue: DRAQ5 fluorescent DNA dye. Red: Actin filaments have been labeled with Alexa Fluor- 555 phalloidin. Secondary antibody from Fisher (Cat#: 35503)




Flow cytometric analysis of HeLa cells using ATPIF1 mouse mAb (green) and negative control (red).

产品名称: ATPIF1 Mouse Monoclonal Antibody
产品货号: AMM81910



Immunohistochemical analysis of paraffin-embedded human lung cancer tissues using ATPIF1 mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human cervical cancer tissues using ATPIF1 mouse mAb with DAB staining.

注意事项 (Note)

For research use only .