

产品名称: HTR3B Mouse Monoclonal Antibody
产品货号: AMM81931



产品概述 (Summary)

产品名称 (Production Name)	HTR3B Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	HTR3B
别名 (Alternative Names)	5-HT3B
基因 ID (Gene ID)	9177.0
蛋白 ID (SwissProt ID)	O95264.Purified recombinant fragment of human HTR3B (AA: extra 23-150) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	50.3kDa

研究背景 (Background)

The product of this gene belongs to the ligand-gated ion channel receptor superfamily. This gene encodes subunit B of the


产品名称: HTR3B Mouse Monoclonal Antibody
产品货号: AMM81931




type 3 receptor for 5-hydroxytryptamine (serotonin), a biogenic hormone that functions as a neurotransmitter, a hormone, and a mitogen. This receptor causes fast, depolarizing responses in neurons after activation. It is not functional as a homomeric complex, but a pentaheteromeric complex with subunit A (HTR3A) displays the full functional features of this receptor.

研究领域 (Research Area)

图片 (Image Data)



Black line: Control Antigen (100 ng); Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng)



Western blot analysis using HTR3B mouse mAb against SH-SY5Y (1) and HeLa (2) cell lysate.

注意事项 (Note)

For research use only .