

产品名称: ADAP1 Mouse Monoclonal Antibody
产品货号: AMM81936



产品概述 (Summary)

产品名称 (Production Name)	ADAP1 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,ELISA,FC
种属反应性 (Reactivity)	Human,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ADAP1
别名 (Alternative Names)	GCS1L; CENTA1; P42ip4
基因 ID (Gene ID)	11033.0
蛋白 ID (SwissProt ID)	O75689 .Purified recombinant fragment of human ADAP1 (AA: 240-370) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	43.4kDa

研究背景 (Background)

ADAP1 (ArfGAP With Dual PH Domains 1) is a Protein Coding gene. Among its related pathways are Class I PI3K signaling


产品名称: ADAP1 Mouse Monoclonal Antibody
产品货号: AMM81936




events and B Cell Receptor Signaling Pathway (sino). GO annotations related to this gene include GTPase activator activity and inositol 1,3,4,5 tetrakisphosphate binding. An important paralog of this gene is ADAP2.

研究领域 (Research Area)


图片 (Image Data)



Black line: Control Antigen (100 ng); Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng)



Western blot analysis using ADAP1 mouse mAb against C6 (1), A549 (2), and HepG2 (3) cell lysate.



Flow cytometric analysis of K562 cells using ADAP1 mouse mAb (green) and negative control (red).

产品名称: ADAP1 Mouse Monoclonal Antibody
产品货号: AMM81936



注意事项 (Note)

For research use only .