

产品名称: P2RY4 Mouse Monoclonal Antibody
产品货号: AMM81956



产品概述 (Summary)

产品名称 (Production Name)	P2RY4 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IHC,ELISA
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	P2RY4
别名 (Alternative Names)	NRU; P2P; UNR; P2Y4
基因 ID (Gene ID)	5030.0
蛋白 ID (SwissProt ID)	P51582 .Purified recombinant fragment of human P2RY4 (AA: extra mix) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:200-1:1000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	41kDa

研究背景 (Background)

The product of this gene belongs to the family of G-protein coupled receptors. This family has several receptor subtypes


产品名称: P2RY4 Mouse Monoclonal Antibody
产品货号: AMM81956




with different pharmacological selectivity, which overlaps in some cases, for various adenosine and uridine nucleotides. This receptor is responsive to uridine nucleotides, partially responsive to ATP, and not responsive to ADP.

研究领域 (Research Area)


图片 (Image Data)



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)




Western blot analysis using P2RY4 mouse mAb against SW620 (1) and A549 (2) cell lysate.



Immunohistochemical analysis of paraffin-embedded human bladder cancer tissues using P2RY4 mouse mAb with DAB staining.

产品名称: P2RY4 Mouse Monoclonal Antibody
产品货号: AMM81956



Immunohistochemical analysis of paraffin-embedded human rectum cancer tissues using P2RY4 mouse mAb with DAB staining.

注意事项 (Note)

For research use only .