

产品名称: CD193 Mouse Monoclonal Antibody
产品货号: AMM81960



产品概述 (Summary)

产品名称 (Production Name)	CD193 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	ELISA,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CD193
别名 (Alternative Names)	CKR3; CCR3; CMKBR3; CC-CKR-3
基因 ID (Gene ID)	1232.0
蛋白 ID (SwissProt ID)	P51677.Purified recombinant fragment of human CD193 (AA: extra mix) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	41kDa

研究背景 (Background)

The protein encoded by this gene is a receptor for C-C type chemokines. It belongs to family 1 of the G protein-coupled


产品名称: CD193 Mouse Monoclonal Antibody
产品货号: AMM81960




receptors. This receptor binds and responds to a variety of chemokines, including eotaxin (CCL11), eotaxin-3 (CCL26), MCP-3 (CCL7), MCP-4 (CCL13), and RANTES (CCL5). It is highly expressed in eosinophils and basophils, and is also detected in TH1 and TH2 cells, as well as in airway epithelial cells. This receptor may contribute to the accumulation and activation of eosinophils and other inflammatory cells in the allergic airway. It is also known to be an entry co-receptor for HIV-1. This gene and seven other chemokine receptor genes form a chemokine receptor gene cluster on the chromosomal region 3p21. Alternatively spliced transcript variants have been described.

研究领域 (Research Area)

图片 (Image Data)




Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Flow cytometric analysis of HL-60 cells using CD193 mouse mAb (green) and negative control (red).

产品名称: CD193 Mouse Monoclonal Antibody
产品货号: AMM81960



Flow cytometric analysis of K562 cells using CD193 mouse mAb (green) and negative control (red).

注意事项 (Note)

For research use only .