

产品名称: CD166 Mouse Monoclonal Antibody
产品货号: AMM81970



产品概述 (Summary)

产品名称 (Production Name)	CD166 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IHC,ELISA,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CD166
别名 (Alternative Names)	MEMD; ALCAM
基因 ID (Gene ID)	214.0
蛋白 ID (SwissProt ID)	Q13740.Purified recombinant fragment of human CD166 (AA: extra 227-381) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:200-1:1000,ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	65kDa

研究背景 (Background)

This gene encodes activated leukocyte cell adhesion molecule (ALCAM), also known as CD166 (cluster of differentiation


产品名称: CD166 Mouse Monoclonal Antibody
产品货号: AMM81970




166), which is a member of a subfamily of immunoglobulin receptors with five immunoglobulin-like domains (VVC2C2C2) in the extracellular domain. This protein binds to T-cell differentiation antigen CD6, and is implicated in the processes of cell adhesion and migration. Multiple alternatively spliced transcript variants encoding different isoforms have been found.

研究领域 (Research Area)


图片 (Image Data)



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)




Flow cytometric analysis of HL-60 cells using CD166 mouse mAb (green) and negative control (red).



产品名称: CD166 Mouse Monoclonal Antibody
产品货号: AMM81970



Immunohistochemical analysis of paraffin-embedded human ovarian cancer tissues using CD166 mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human rectum cancer tissues using CD166 mouse mAb with DAB staining.

注意事项 (Note)

For research use only .