

产品名称: P2RY12 Mouse Monoclonal Antibody
产品货号: AMM81977



产品概述 (Summary)

产品名称 (Production Name)	P2RY12 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,ELISA,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	P2RY12
别名 (Alternative Names)	HORK3; P2Y12; ADPG-R; BDPLT8; SP1999; P2T(AC); P2Y(AC); P2Y(12)R; P2Y(ADP); P2Y(cyc)
基因 ID (Gene ID)	64805.0
蛋白 ID (SwissProt ID)	Q9H244.Purified recombinant fragment of human P2RY12 (AA: extra mix) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	39.4kDa

研究背景 (Background)


产品名称: P2RY12 Mouse Monoclonal Antibody
产品货号: AMM81977




The product of this gene belongs to the family of G-protein coupled receptors. This family has several receptor subtypes with different pharmacological selectivity, which overlaps in some cases, for various adenosine and uridine nucleotides. This receptor is involved in platelet aggregation, and is a potential target for the treatment of thromboembolisms and other clotting disorders. Mutations in this gene are implicated in bleeding disorder, platelet type 8 (BDPLT8). Alternative splicing results in multiple transcript variants of this gene.

研究领域 (Research Area)

图片 (Image Data)




Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Western blot analysis using P2RY12 mouse mAb against PC-3 (1) and C6 (2) cell lysate.

产品名称: P2RY12 Mouse Monoclonal Antibody
产品货号: AMM81977



Flow cytometric analysis of Hela cells using P2RY12 mouse mAb (green) and negative control (red).

注意事项 (Note)

For research use only .