

产品名称: PLXNA1 Mouse Monoclonal Antibody
产品货号: AMM81981



产品概述 (Summary)

产品名称 (Production Name)	PLXNA1 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	ICC,ELISA,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PLXNA1
别名 (Alternative Names)	NOV; NOVP; PLXN1; PLEXIN-A1
基因 ID (Gene ID)	5361.0
蛋白 ID (SwissProt ID)	Q9UIW2.Purified recombinant fragment of human PLXNA1 (AA: 1100-1200) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	ICC 1:100-1:500,ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	211kDa

研究背景 (Background)

PLXNA1 (Plexin A1) is a Protein Coding gene. Diseases associated with PLXNA1 include Hereditary Congenital Facial Paresis.


产品名称: PLXNA1 Mouse Monoclonal Antibody
产品货号: AMM81981




Among its related pathways are Developmental Biology and Neuroscience. GO annotations related to this gene include receptor activity and semaphorin receptor activity. An important paralog of this gene is PLXNA3.

研究领域 (Research Area)


图片 (Image Data)



Black line: Control Antigen (100 ng); Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng)



Immunofluorescence analysis of HeLa cells using PLXNA1 mouse mAb (green). Blue: DRAQ5 fluorescent DNA dye. Red: Actin filaments have been labeled with Alexa Fluor- 555 phalloidin. Secondary antibody from Fisher (Cat#: 35503)



Flow cytometric analysis of HeLa cells using PLXNA1 mouse mAb (green) and negative control (red).

产品名称: PLXNA1 Mouse Monoclonal Antibody
产品货号: AMM81981



注意事项 (Note)

For research use only .