

产品名称: TIM3 Mouse Monoclonal Antibody
产品货号: AMM81997



产品概述 (Summary)

产品名称 (Production Name)	TIM3 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	ICC,ELISA,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TIM3
别名 (Alternative Names)	HAVCR2; CD366; KIM-3; TIMD3; Tim-3; TIMD-3; HAVcr-2
基因 ID (Gene ID)	84868.0
蛋白 ID (SwissProt ID)	Q8TDQ0.Purified recombinant fragment of human TIM3 (AA: extra 22-202) expressed in HEK293 cells.

产品应用 (Application)

稀释比 (Dilution Ratio)	ICC 1:200-1:1000,ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	33.4kDa

研究背景 (Background)

The protein encoded by this gene belongs to the immunoglobulin superfamily, and TIM family of proteins. CD4-positive T


产品名称: TIM3 Mouse Monoclonal Antibody
产品货号: AMM81997




helper lymphocytes can be divided into types 1 (Th1) and 2 (Th2) on the basis of their cytokine secretion patterns. Th1 cells are involved in cell-mediated immunity to intracellular pathogens and delayed-type hypersensitivity reactions, whereas, Th2 cells are involved in the control of extracellular helminthic infections and the promotion of atopic and allergic diseases. This protein is a Th1-specific cell surface protein that regulates macrophage activation, and inhibits Th1-mediated auto- and alloimmune responses, and promotes immunological tolerance.

研究领域 (Research Area)

图片 (Image Data)




Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Immunofluorescence analysis of HeLa cells using TIM3 mouse mAb (green). Blue: DRAQ5 fluorescent DNA dye. Red: Actin filaments have been labeled with Alexa Fluor- 555 phalloidin. Secondary antibody from Fisher (Cat#: 35503)

产品名称: TIM3 Mouse Monoclonal Antibody
产品货号: AMM81997



Flow cytometric analysis of Jurkat cells using TIM3 mouse mAb (green) and negative control (red).

注意事项 (Note)

For research use only .