

产品名称: VISTA Mouse Monoclonal Antibody
产品货号: AMM81998



产品概述 (Summary)

产品名称 (Production Name)	VISTA Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	ELISA,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	VISTA
别名 (Alternative Names)	VSIR; B7H5; GI24; B7-H5; PD-1H; SISP1; PP2135; C10orf54; DD1alpha
基因 ID (Gene ID)	64115.0
蛋白 ID (SwissProt ID)	Q9H7M9.Purified recombinant fragment of human VISTA (AA: extra 33-194) expressed in HEK293 cells.

产品应用 (Application)

稀释比 (Dilution Ratio)	ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	34kDa

研究背景 (Background)

VSIR (V-Set Immunoregulatory Receptor) is a Protein Coding gene. Diseases associated with VSIR include Parasagittal


产品名称: VISTA Mouse Monoclonal Antibody
产品货号: AMM81998




Meningioma and Solar Retinopathy. Among its related pathways are T Cell Co-Signaling Pathway: Ligand-Receptor Interactions and NF-kappaB Signaling.

研究领域 (Research Area)

图片 (Image Data)



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Flow cytometric analysis of Jurkat cells using VISTA mouse mAb (green) and negative control (red).

注意事项 (Note)

For research use only .