

产品名称: APAF1 Mouse Monoclonal Antibody
产品货号: AMM82023



产品概述 (Summary)

产品名称 (Production Name)	APAF1 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	ELISA,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG2b
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	APAF1
别名 (Alternative Names)	CED4; APAF-1
基因 ID (Gene ID)	317.0
蛋白 ID (SwissProt ID)	O14727 .Purified recombinant fragment of human APAF1 (AA: 1138-1237) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	141.8kDa

研究背景 (Background)

This gene encodes a cytoplasmic protein that initiates apoptosis. This protein contains several copies of the WD-40 domain,


产品名称: APAF1 Mouse Monoclonal Antibody
产品货号: AMM82023




a caspase recruitment domain (CARD), and an ATPase domain (NB-ARC). Upon binding cytochrome c and dATP, this protein forms an oligomeric apoptosome. The apoptosome binds and cleaves caspase 9 preproprotein, releasing its mature, activated form. Activated caspase 9 stimulates the subsequent caspase cascade that commits the cell to apoptosis. Alternative splicing results in several transcript variants encoding different isoforms.

研究领域 (Research Area)

图片 (Image Data)



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Flow cytometric analysis of HeLa cells using APAF1 mouse mAb (green) and negative control (red).

注意事项 (Note)

For research use only .