

产品名称: CD11C Mouse Monoclonal Antibody
产品货号: AMM82034



产品概述 (Summary)

产品名称 (Production Name)	CD11C Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	ELISA,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CD11C
别名 (Alternative Names)	ITGAX; SLEB6
基因 ID (Gene ID)	3687.0
蛋白 ID (SwissProt ID)	P20702.Purified recombinant fragment of human CD11C (AA: extra 102-279) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	127.8kDa

研究背景 (Background)

This gene encodes the integrin alpha X chain protein. Integrins are heterodimeric integral membrane proteins composed of


产品名称: CD11C Mouse Monoclonal Antibody
产品货号: AMM82034




an alpha chain and a beta chain. This protein combines with the beta 2 chain (ITGB2) to form a leukocyte-specific integrin referred to as inactivated-C3b (iC3b) receptor 4 (CR4). The alpha X beta 2 complex seems to overlap the properties of the alpha M beta 2 integrin in the adherence of neutrophils and monocytes to stimulated endothelium cells, and in the phagocytosis of complement coated particles. Two transcript variants encoding different isoforms have been found for this gene.

研究领域 (Research Area)

图片 (Image Data)



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Flow cytometric analysis of Raji cells using CD11C mouse mAb (green) and negative control (red).

注意事项 (Note)

For research use only .