

产品名称: CD84 Mouse Monoclonal Antibody
产品货号: AMM82052



产品概述 (Summary)

产品名称 (Production Name)	CD84 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	ELISA,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CD84
别名 (Alternative Names)	LY9B; hCD84; mCD84; SLAMF5
基因 ID (Gene ID)	8832.0
蛋白 ID (SwissProt ID)	Q9UIB8.Purified recombinant fragment of human CD84 (AA: extra 22-225) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	38.8kDa

研究背景 (Background)

This gene encodes a membrane glycoprotein that is a member of the signaling lymphocyte activation molecule (SLAM)

产品名称: CD84 Mouse Monoclonal Antibody
产品货号: AMM82052




family. This family forms a subset of the larger CD2 cell-surface receptor Ig superfamily. The encoded protein is a homophilic adhesion molecule that is expressed in numerous immune cells types and is involved in regulating receptor-mediated signaling in those cells. Alternate splicing results in multiple transcript variants.


研究领域 (Research Area)

Autophagy

图片 (Image Data)



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Flow cytometric analysis of HL-60 cells using CD84 mouse mAb (green) and negative control (red).

注意事项 (Note)

For research use only .