

产品名称: TLL1 Mouse Monoclonal Antibody
产品货号: AMM82057



产品概述 (Summary)

产品名称 (Production Name)	TLL1 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	ELISA,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TLL1
别名 (Alternative Names)	TLL; ASD6
基因 ID (Gene ID)	7092.0
蛋白 ID (SwissProt ID)	O43897.Purified recombinant fragment of human TLL1 (AA: 870-1013) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	114.7kDa

研究背景 (Background)

This gene encodes an astacin-like, zinc-dependent, metalloprotease that belongs to the peptidase M12A family. This


产品名称: TLL1 Mouse Monoclonal Antibody
产品货号: AMM82057




protease processes procollagen C-propeptides, such as chordin, pro-biglycan and pro-lysyl oxidase. Studies in mice suggest that this gene plays multiple roles in the development of mammalian heart, and is essential for the formation of the interventricular septum. Allelic variants of this gene are associated with atrial septal defect type 6. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.

研究领域 (Research Area)

图片 (Image Data)



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Flow cytometric analysis of HL-60 cells using TLL1 mouse mAb (green) and negative control (red).

注意事项 (Note)

For research use only .